

EverLift[®]
Submucosal Lifting Agent




CONVENIENTLY
Packaged


SUPERIOR
Cushion


FOR EVERYDAY
Use



NO PREPARATION REQUIRED
AS EASY TO USE AS SPOT[®] EX ENDOSCOPIC TATTOO



Pre-mixed with methylene blue
→ saving time vs hospital-mixed solutions



Supplied **sterile** in a pre-filled syringe
→ reducing infection risk



Stored at **room temperature**
→ no refrigeration is required




Syringe attaches directly to **any endoscopic injection needle's luer-lock connection** (22G/23G needle recommended)



Conforms to MHR guidelines¹ on the "restricted use of injectable medication"
→ minimising risk to patients, staff & supply chain

1. NHSI Patient Safety Alert NHS/PSA/D/2016/008

INDICATION FOR USE: EverLift Submucosal Lifting Agent is indicated for use in gastrointestinal endoscopic procedures for submucosal lift of polyps, adenomas, early stage cancers, or other gastrointestinal lesions prior to excision with a snare or other appropriate endoscopic device



ENDOMED PTY. LTD.
SPECIALISING IN ENDOSCOPY

A.B.N. 98 050 413 379

orders@endomed.com.au

ENDOMED PTY. LTD.
POST OFFICE BOX 445,
STRATHPINE QLD. 4500

TELEPHONE 07 3881 1883
1800 773 474
FACSIMILE 07 3881 0211


THE LIFT you've been waiting for



SUPERIOR CUSHION



FOR EVERYDAY USE



Cellulose based & Hypertonic: with a slower rate of absorption due to osmotic resistance

EverLift provides margin delineation with darker color



EverLift sustains a 2mm lift for 60 minutes²

2. Data on file

Why choose EverLift?

"From a nursing perspective, EverLift is an ideal and very efficient lifting agent. It is not as viscous as some of the other agents, so it is very easy for the nurse to inject. During difficult endoscopic mucosal resections, EverLift keeps the submucosal layer lifted for longer which requires fewer injections. This ensures a more seamless resection because it requires fewer steps. It is also more cost-effective than other lifting agents. It requires less volume to lift effectively, and the 5 mL syringes are an excellent way to eliminate waste which makes it a great option for any hospital. EverLift is a great tool to have in your arsenal!"

Would you recommend EverLift to a colleague?

"I would recommend EverLift to a colleague without reservation. The nurses and physicians at AUMC have been very pleased with this product!"

Bonnie Clark, RN

Advanced Endoscopy Registered Nurse
Augusta University Medical Center-Digestive Health Center, Augusta



ORDER NO.	DESCRIPTION	QUANTITY
GIS-55	EverLift Submucosal Lifting Agent 5 mL	Box of 10 Syringes
GIS-59	EverLift Submucosal Lifting Agent 10 mL	Box of 10 Syringes

Disclaimer: EverLift Submucosal Lifting Agent is FDA approved and CE Marked. Please check availability with your local GI Supply distributor.

ENDOMED PTY. LTD.
SPECIALISING IN ENDOSCOPY

A.B.N. 98 050 413 379

orders@endomed.com.au

ENDOMED PTY. LTD.
POST OFFICE BOX 445,
STRATHPINE QLD. 4500

TELEPHONE 07 3881 1883
1800 773 474
FACSIMILE 07 3881 0211