

# Bite Blocks

AUSTRALIAN DESIGNED RANGE OF DISPOSABLE BITE BLOCKS

VARIETY OF STYLES TO SUIT YOUR PATIENT'S NEEDS



DISTRIBUTED EXCLUSIVELY THROUGHOUT AUSTRALIA BY

**ENDOMED PTY. LTD.**

**1800 773 474**

07 3881 1883

MEDIKAL PRODUCTS' EXTENSIVE RANGE OF BITE BLOCKS OFFER A VARIETY OF MODELS TO COVER THE NEEDS OF MODERN ENDOSCOPY, ENSURING THE BEST OUTCOME FOR PATIENT AND PRACTITIONER

## KEY FEATURES AND BENEFITS

### STREAMLINED DESIGN – MAXIMUM COMFORT

Different models for clinical versatility

### RIGID INNER CORE

Protecting your device from potential “bite-through” effect

### READY TO USE

Each Bite Block is individually wrapped in a quick-to-open bag

### SINGLE USE AND LATEX FREE

No expensive reprocessing

### TRACKING STICKERS

Individual Patient Tracking Stickers included



### GlideGuard

KM-0201

- 47 Fr Disposable Bite Block
- Sculpted Padding for patients with compromised mouths
- Smooth Strapless Design



### GlideGuard

KM-0202

- 47 Fr Disposable Bite Block
- Sculpted Padding for patients with compromised mouths
- Velcro Head Strap included



### EasyGuard

KM-0203

- 48 Fr Disposable Bite Block
- Streamlined design
- Latex Free Head Strap included
- Dispensing carton



### MaxGuard

KM-0204

- 72 Fr Disposable Bite Block
- Wide flanges for optimal positioning
- Latex Free Head Strap included
- Dental retention ring
- Dispensing carton



### PediGuard

KM-0205

- 39 Fr Disposable Bite Block
- Streamlined design for infants and smaller mouths
- Latex Free Head Strap included



### GenGuard

KM-0206

- 60 Fr Disposable Bite Block
- Smooth Strapless Design
- Dental retention ring
- Dispensing carton



### PediGuard Soft

KM-0215

- 39 Fr Disposable Padded Bite Block
- Sculpted Padding for patients with compromised mouths
- Streamlined design for infants and smaller mouths
- Latex Free Head Strap included

PRODUCT NUMBER	SIZE (INTERNAL OPENING)	DESCRIPTION	UNIT OF MEASURE
KM-0201	47 Fr (15.7mm x 28mm)	GlideGuard (Strapless)	Box of 24
KM-0202	47 Fr (15.7mm x 28mm)	GlideGuard With Velcro Strap	Box of 24
KM-0203	48 Fr (16mm x 26mm)	EasyGuard with Latex Free Strap	Box of 100
KM-0204	72 Fr (24mm x 27mm)	MaxGuard with Latex Free Strap	Box of 100
KM-0205	39 Fr (13mm x 20mm)	PediGuard with Latex Free Strap	Box of 25
KM-0215	39 Fr (13mm x 20mm)	PediGuard Soft with Latex Free Strap	Box of 25
KM-0206	60 Fr (20mm x 28mm)	GenGuard (Strapless)	Box of 100

For more information or samples please contact

**ENDOMED PTY. LTD.** 1800 773 474  
info@endomed.com.au 07 3881 1883